MAS® UA Control

LIQUID ASSAYED URINALYSIS CONTROL

IVD

INTENDED USE

MAS® UA Control is intended for use in the clinical laboratory as a control for qualitative and semi-quantitative procedures used in routine urinalysis testing. Assay values are provided for specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and equipment. The product is suitable for use as a control material for physiochemical, chemical, and microscopic methods of routine urinalysis. **UA Control** may be used in conjunction with commercially available standardized systems for urine microscopic analysis.

PRODUCT DESCRIPTION

UA Control is a liquid stable control material prepared from human urine. Analyte levels are adjusted with various pure chemicals and human source materials. **UA Control** also contains preservatives and stabilizers.

UA Control is prepared from human urine pools and contains other human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antibody (HCV) and HIV-1 and HIV-2 Antibody. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from the Center for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 2009.

DANGER: UA Control contains <0.9% bovine serum albumin (BSA), <0.1% Sulfobromophthalein sodium salt and <0.006% 3:1 Mixture: 5-chloro-2-methyl-4-isothiazolin-3-one + 2-methyl-4-isothiazolin-3-one (PROCLIN)

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

STORAGE AND STABILITY

Store unopened **UA Control** at 2-8°C. Unopened vials are stable until the expiration date on the label. Once opened, vials of control are stable for 6 weeks when stored tightly capped at room temperature (18-25°C) or 3 months refrigerated (2-8°C). **Do not freeze.**

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

EXPECTED VALUES

Expected values may vary slightly with different lot numbers and reagent strips. Refer to the included table for values obtained on specific systems. Good laboratory practice suggests that each laboratory establish its own parameters.

Values listed were obtained using instrument(s) and instrument manufacturer's reagents and procedures at the time of assay. Subsequent modifications in instrument, reagent, or procedures may invalidate these values.

Laboratories using methods other than those listed should establish their own expected ranges.

INSTRUCTIONS FOR USE

Materials Required But NOT Provided:

- Slides
- Pipettors
- Centrifuge tubes (Centrifuge tubes, slides, and pipettors from standardized urinalysis systems may be used.)
- Centrifuge capable of generating an RCF of 400 for KOVA method or 450 for manual methods.
- Microscope with a 10x objective and a 40x objective lens
- Urinalysis stain (Optional)
- Chemical reagent strips
- Any additional laboratory equipment or reagents required to perform the physiochemical, chemical and microscopic analyses required.

Allow the **UA Control** to come to room temperature. Thoroughly mix the contents of the vial before each use by gently inverting for several minutes (3 minutes minimum). Replace cap immediately and store the opened vial at room temperature.

Physical Tests:

1. **Appearance:** Visually inspect the color and turbidity.
2. **Specific Gravity:** Measure the specific gravity using a refractometer or reagent strip.
3. **Osmolality:** Measure osmolality with an osmometer.

Chemical Tests:

1. If reading results visually, compare the reagent strip with the manufacturer's color comparison chart. If using an automated strip reader, follow manufacturer's directions for use.

2a. **15 mL Bottle:** Remove cap from bottle and pour an aliquot of control into a clean test tube. Replace cap immediately after use. Immerse the reagent strip into the control. Remove the strip and follow manufacturer's instructions for testing.

2b. **60 mL Bottle:** Open the vial and pour an aliquot of control into a clean test tube. Replace cap immediately after use. Immerse the reagent strip into the control. Remove the strip and follow manufacturer's instructions for testing.

Microscopic Analysis:

If using a standardized methodology for microscopic analysis, follow manufacturer's recommendations. Each laboratory should establish its own values for systems other than those listed.

If using a non-standardized methodology for microscopic analysis, use a 12 mL sample centrifuged for 5 minutes at 450 RCF (relative centrifugal force). Pour off supernatant leaving approximately 0.5 mL. Gently resuspend sediment transferring one drop onto a clean, dry microscope slide and cover with a 22 mm coverslip. Values reported are the range of microscopic elements counted in 10 hpf (high powered fields). Each laboratory should establish its own ranges.

Handle the **UA Control** in the same manner as patient specimens. Consistent handling of the control is necessary for uniform results.

This product contains cellular elements which may form a small amount of sediment in the bottle. To ensure consistent recoveries, invert and mix well prior to use.

15 mL Bottle: Wearing eye protection, carefully remove and discard the dropper tip from the bottle. Pour 12 mL of control into a clean conical centrifuge tube. A consistent volume for testing is essential for reproducible microscopic results. Store opened vial at room temperature.

60 mL Bottle: Open the vial and pour 12 mL of control into a clean conical centrifuge tube. A consistent volume for testing is essential for reproducible microscopic results. Replace cap immediately after use. Store opened vial at room temperature.

QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURE

UA Control must be mixed well before each use to ensure reproducible results.

A consistent volume must be used for testing to ensure reproducible microscopic results.

Organic spheres used to simulate leukocytes may not exhibit the same staining characteristics of human leukocytes.

See manufacturer's product insert for a listing of limitations for each test method. For professional use only.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit www.maslablink.com and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit www.maslablink.com select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
UAB-115	UA Control, Normal, Level 1	6 x 15 mL
UAB-215	UA Control, Abnormal, Level 2	6 x 15 mL
UAB-MP	UA Control, Multi-Pack	6 x 15 mL
UAB-160	UA Control, Normal, Level 1	4 x 60 mL
UAB-260	UA Control, Abnormal, Level 2	4 x 60 mL

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LVL 1 2026-01-31
LVL 2 2026-01-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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

For Insert updates go to:
www.thermofisher.com/diagnostics

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Acetest®, Advantus, Atlas, Clinitek® Series, Clinitest®, Ictotest®, Micro-Bumintest®, Status®, Diastix®, Reg. TM: Siemens Healthcare Diagnostic, Tarrytown, NY
Array 360®, ICON®, Image®, Reg. TM: Beckman Coulter Inc., Fullerton, CA
Biovation, Reg. TM: Biovation, Inc., Hercules, CA
Dimension®, Reg. TM: Siemens Healthcare Diagnostic, Glasgow, DE
URI-system®, Sure-View, Reg. TM: Thermo Fisher Scientific, Waltham, MA
Signify®, Contrast, OSOM Reg. TM: Genzyme Corporation, Cambridge, MA
KOVA®, Reg. TM: Hycor Biomedical Inc., Garden Grove, CA
DiaScreen Series, Reg. TM: Hypoguard, Minneapolis, MN
Iris IQ200, Yellow IRIS®, Reg. TM: International Remote Imaging Systems, Inc., Chatsworth, CA
Chemstrip®, Chemstrip® 101, Criterion, CUA, Hitachi®, Micral, Mini UA, SUA, URISYS 1100, URISYS Systems, Reg. TM: Roche Diagnostics, Indianapolis, IN
QuickVue®, UrinQuick, Reg. TM: Quidel Corporation, San Diego, CA
Cen-Slide®, Reg. TM: StatSpin, Norwood, MA
PocketChem UA®, Reg. TM: ThermoFisher Scientific, Boulder, CO
YD Uriscan®, Reg. TM: YD Diagnostics, Seoul, Korea

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


If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.

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USA
☎ 800-232-3342
sales.diagnostics.fmt@thermofisher.com

China
☎ +86 800-810-5118
cdx.cn.info@thermofisher.com

New Zealand
☎ 0800 933 966
auinfo@thermofisher.com

Switzerland & Austria
☎ +41 26 663 86 70
cdx.ch.info@thermofisher.com

Asia Pacific
☎ +61 1800 333 110
cdd.asia.info@thermofisher.com

France
☎ +33 1 40 86 65 20
cdx.fr.info@thermofisher.com

Nordic
☎ +47 2 325 0433
info.nordic.cdd@thermofisher.com

United Kingdom & Ireland
☎ +44 1442 868 940
cdx.uk.info@thermofisher.com

Australia
☎ +61 1800 333 110
auinfo@thermofisher.com

Germany
☎ +49 0800-40 40 771
cdx.de.info@thermofisher.com

South Africa
☎ +27117926790
support-za.idd@thermofisher.com

Canada
☎ 800-282-4075
info.cddcanada@thermofisher.com

Japan
☎ +81 (0)120-147-075
JPYOK-CDD.QC@thermofisher.com

Spain, Portugal & Italy
☎ +34 93589 8338
cdx.es.info@thermofisher.com

For countries not listed:
distributor.cdd@thermofisher.com

Microgenics Corporation
46500 Kato Road
Fremont, CA 94538 USA
ISO 13485 Certified Company
www.thermofisher.com/diagnostics

EC REP

B.R.A.H.M.S GmbH
Neuendorffstrasse 25
16761 Hennigsdorf, Germany

UAB-INS-VA
Rev. 22 2020 11

thermo
scientific