MAS® Urinalysis

LVL 1

LIQUID ASSAYED URINALYSIS CONTROL

IVD

INTENDED USE

MAS® Urinalysis is intended for use in the clinical laboratory as a control for qualitative and semi-quantitative procedures used in routine urinalysis testing. Assay values are provided for specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and equipment. The product is suitable for use as a control material for physiochemical, chemical, and microscopic methods of routine urinalysis. **Urinalysis** may be used in conjunction with commercially available standardized systems for urine microscopic analysis.

PRODUCT DESCRIPTION

Urinalysis is a liquid stable control material prepared from human urine. Analyte levels are adjusted with various pure chemicals and preparations of human erythrocytes, human leukocytes and hCG from human source material, and simulated leukocytes. **Urinalysis** also contains preservatives and stabilizers.

NOTE: MICROALBUMIN TESTING

Urinalysis is formulated for use in microalbumin testing. To do so, the user should assay **Urinalysis L-3** as a normal negative. Although **Urinalysis L-1** contains albumin, this concentration is relatively high compared to abnormal ranges for microalbumin. Consequently, results with most reagent manufacturer's products will yield the highest reportable level.

CAUTION: Urinalysis is prepared from human urine pools and contains other human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antibody (HCV) and HIV-1 and HIV-2 Antibody. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from the Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 2009. The packaging of this product contains dry natural rubber.

DANGER: Urinalysis control contains ≤0.67% Leukocyte Esterase solution (LE) and ≤0.39% human serum albumin (HSA).

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/ face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

STORAGE AND STABILITY

Store **Urinalysis** at 2-8°C. Unopened vials are stable until the expiration date on the label. Once opened, vials of control are stable for 30 days when stored tightly capped at 2-8°C. **Do not freeze.**

Urobilinogen (**Urinalysis, Level 1**) is subject to degradation when exposed to air. Use as directed to ensure product stability.

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

CONTROL RANGES

The published control ranges are based upon a combination of replicate assays of representative samples by participating laboratories, instrument/reagent manufacturers and direct correlation with other analytical systems in accordance with established protocol. Instrument values provided are specific to this lot of control only and are intended to assist the laboratory in establishing its own means and ranges. All values have been assigned with instruments and reagents available at the time of assay and expected values may vary with different reagents and/or methodologies. Laboratory established means should fall within the assigned ranges although subsequent instrument, reagent or calibration modifications may invalidate assigned values.

Peer comparison data and latest QC lot specific updates are available online through LabLink® xL Quality Assurance Program at www.maslablink.com. Refer to the Technical Assistance section for contact information.

INSTRUCTIONS FOR USE

Materials Required But NOT Provided:

- Slides
- Pipettors
- Centrifuge tubes (Centrifuge tubes, slides, and pipettors from standardized urinalysis systems may be used.)

- Centrifuge capable of generating an RCF of 400 for Kova method or 450 for manual methods.
- Microscope with a 10X objective and a 40X objective lens
- Urinalysis stain (Optional)
- Chemical reagent strips
- Any additional laboratory equipment or reagents required to perform the physiochemical, chemical and microscopic analyses required.

Thoroughly mix the contents of the vial before each use by gently inverting for several minutes (3 minutes minimum). Open the vial and pour 12.0 mL of control into a clean conical centrifuge tube. Replace cap immediately and store the opened vial at 2-8°C. A consistent volume for testing is essential for reproducible microscopic results.

Allow the control in the tube to come to room temperature (20-25°C) before testing.

Physical Tests:

- Appearance:** Record the color and turbidity in the appropriate column on a daily Quality Control Chart.
- Specific Gravity:** Measure the specific gravity using a refractometer, urinometer, or reagent strip and record results on a daily Quality Control Chart.
- Osmolality:** Measure osmolality with an osmometer. Record results in the appropriate column.

Chemical Tests:

- Follow the manufacturer's directions for the specific chemistry reagent strip.
- Immerse the reagent strip into the urine. Remove the strip and blot to remove the excess urine.
- Visually compare the reagent strip with the manufacturer's color comparison chart or read using an automated strip reader following manufacturer's directions for use.

Microscope Analysis:

If using a standardized methodology for microscopic analysis, follow manufacturer's recommendations. Each laboratory should establish its own values for systems other than those listed.

If using a non-standardized methodology for microscopic analysis, use a 12 mL sample centrifuged for 5 minutes at 450 RCF (relative centrifugal force). Pour off supernatant leaving approximately 0.5 mL. Gently resuspend sediment transferring one drop onto a clean, dry microscope slide and cover with a 22 mm coverslip. Scan the slide chamber to enumerate microscopic elements in 10 hpf (high powered fields). Each laboratory should establish its own ranges.

Handle the **Urinalysis** in the same manner as patient specimens. Consistent handling of the control is necessary for uniform results.

This product contains cellular elements which may form a small amount of sediment in the bottle. To ensure consistent recoveries, invert and mix well prior to use.

QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURE

Urinalysis must be mixed well before each use to ensure reproducible results.

A consistent volume must be used for testing to ensure reproducible microscopic results.

Organic spheres used to simulate leukocytes may not exhibit the same staining characteristics of human leukocytes.

See manufacturer's product insert for a listing of limitations for each test method. For professional use only.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit www.maslablink.com and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit www.maslablink.com select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
024222	Urinalysis, High abnormal, Level 1	4 x 15 mL
024224	Urinalysis, Normal, Level 3	4 x 15 mL
024225	Urinalysis, High abnormal, Level 1	4 x 60 mL
024227	Urinalysis, Normal, Level 3	4 x 60 mL

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
IVD

		 UA1525031/UA6025031
Roche Urisys 2400	Specific Gravity pH Leukocyte Esterase Nitrite Protein Glucose Ketones Urobilinogen Bilirubin Blood (Hemoglobin)	1.005 - 1.050 7.0 - 9.0 25(+) - 500(+++) Leuk/μL Positive 25(+) - 500(+++++) mg/dL or 0.25 - 5.00 g/L 100(++) - 1000(+++++) mg/dL or 6 - 56 mmol/L 5(+) - 150(+++++) mg/dL or 0.5 - 15 mmol/L 1(+) 12(+++++) mg/dL (EU/dL) or 17 - 203 μmol/L 1(+) - 6(+++++) mg/dL or 17 - 100 μmol/L 10(+) - 250(+++++) Ery/μL
Siemens Reagent Strip (Visual)	Leukocyte Esterase Nitrite Urobilinogen Protein pH Blood (Hemoglobin) Specific Gravity Ketones Bilirubin Glucose Creatinine Protein-to-Creatinine ratio	Trace - Large Positive 2.0 - 8.0 EU/dL 100 - 2000 mg/dL 7.0 - 8.5 Moderate - Large 1.005 - 1.025 5 - 160 mg/dL Small - Large 100 - 2000 mg/dL 100 - 300 mg/dL Abnormal
Siemens Clinitek[®] 50	Glucose Bilirubin Ketones Specific Gravity Blood (Hemoglobin) pH Protein Urobilinogen Nitrite Leukocyte Esterase Microalbumin Creatinine Protein-to-Creatinine ratio	100 - ≥ 1000 mg/dL Small - Large 5 - ≥ 80 mg/dL 1.005 - 1.025 Moderate - Large 7.0 - ≥ 9.0 100 - ≥ 300 mg/dL 2.0 - ≥ 8.0 EU/dL Positive Trace - Large ≤ 150 mg/dL 100 - 300 mg/dL Abnormal
Siemens Clinitek[®] 100	Glucose Bilirubin Ketones Specific Gravity Blood (Hemoglobin) pH Protein Urobilinogen Nitrite Leukocyte Esterase Microalbumin Creatinine Protein-to-Creatinine ratio	Trace - ≥ 1000 mg/dL Small - Large 5 - ≥ 80 mg/dL ≤1.005 - 1.025 Moderate - Large 7.0 - ≥ 9.0 100 - ≥ 300 mg/dL 2.0 - ≥ 8.0 EU/dL Positive Trace - Large ≥ 150 mg/L 100 - 300 mg/dL Abnormal
Siemens Clinitek[®] 500/ Clinitek[®] Advantus	Glucose Bilirubin Ketones Specific Gravity Blood (Hemoglobin) pH Protein Urobilinogen Nitrite Leukocyte Esterase Creatinine Protein-to-Creatinine ratio	100(Trace) - ≥ 1000 mg/dL or 6 - ≥ 56 mmol/L Small(+) - Large(+++) Trace(Trace) - ≥ 80(+++++) mg/dL or 0.5 - ≥ 8.0 mmol/L ≤1.005 - 1.025 Small(+) - Large(+++) or Ca 25 - Ca 200 Ery/μL 7.0 - ≥ 9.0 100(++) - ≥300(+++++) mg/dL or 1 - ≥3 g/L 1.0 - ≥ 8.0 EU/dL Positive Trace(Trace) - Large(+++) or Ca 15 - Ca 500 Leuk/μL 100 - 300 mg/dL or 8.84 - 26.5 mmol/L Abnormal
Siemens Clinitek[®] Atlas	Glucose Bilirubin Ketones Specific Gravity Blood (Hemoglobin) pH Protein Urobilinogen Nitrite Leukocyte Esterase	100(Trace) - ≥1000(+++++) mg/dL or 5.6 - ≥56 mmol/L Small - Large 5 - ≥ 80 mg/dL 1.020 - ≥ 1.035 Small - Large 7.0 - 8.5 100 - ≥ 300 mg/dL 1.0 - ≥ 8.0 EU/dL Positive Small - Large

LOT UA1525031
UA6025031

 2025-03-31

thermo
scientific


For insert updates go to:
www.thermofisher.com/diagnostics

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

If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.

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USA
☎ 800-232-3342
sales.diagnostics.fmt@thermofisher.com

China
☎ +86 800-810-5118
cdx.cn.info@thermofisher.com

New Zealand
☎ 0800 933 966
auinfo@thermofisher.com

Switzerland & Austria
☎ +41 26 663 86 70
cdx.ch.info@thermofisher.com

Asia Pacific
☎ +61 1800 333 110
cdd.asia.info@thermofisher.com

France
☎ +33 1 40 86 65 20
cdx.fr.info@thermofisher.com

Nordic
☎ +47 2 325 0433
info.nordic.cdd@thermofisher.com

United Kingdom & Ireland
☎ +44 1442 868 940
cdx.uk.info@thermofisher.com

Australia
☎ +61 1800 333 110
auinfo@thermofisher.com

Germany
☎ +49 0800-40 40 771
cdx.de.info@thermofisher.com

South Africa
☎ +27117926790
support-za.idd@thermofisher.com

Canada
☎ 800-282-4075
info.cddcanada@thermofisher.com

Japan
☎ +81 (0)120-147-075
JPYOK-CDD.QC@thermofisher.com

Spain, Portugal & Italy
☎ +34 93589 8338
cdx.es.info@thermofisher.com

For countries not listed:
distributor.cdd@thermofisher.com

Microgenics Corporation
46500 Kato Road
Fremont, CA 94538 USA
ISO 13485 Certified Company
www.thermofisher.com/diagnostics

B.R.A.H.M.S GmbH
Neuendorfstrasse 25
16761 Hennigsdorf, Germany
Tel: +49 (0) 800 404 077 12
Fax: +49 (0) 800 404 077 13

UA1-INS-VA
Rev. 56 2020 12

thermo
scientific