MAS® Omni•CARDIO™

LIQUID ASSAYED INTEGRATED CARDIAC CONTROL

IVD

Read Highlighted Changes: Revised March 2023

INTENDED USE

Thermo Scientific MAS® Omni•CARDIO™ is intended for use in the clinical laboratory as an assayed control serum for monitoring assay conditions related to cardiac and associated critical marker determinations.

PRODUCT DESCRIPTION

Omni•CARDIO is a liquid stable control material prepared from human serum. Analyte levels are adjusted with various pure chemicals and preparations from recombinant proteins, human tissue or body fluids. Preservatives and stabilizers are added to maintain product integrity.

CAUTION: Omni•CARDIO is prepared from human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C (HCV), HIV-1 and HIV-2. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories", 2009. The packaging of this product contains dry natural rubber.

DANGER: Omni•CARDIO control contains ≤4.0% Human Serum Albumin (HSA), ≤0.02% Chlorhexidine diacetate and ≤0.02% Prionex Porcine Gel.

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H412 - Harmful to aquatic life with long lasting effects.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Avoid release to the environment. Wear protective gloves/eye protection/face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

CONSTITUENT LIST

The following constituents are contained in Omni•CARDIO:

Beta-Human Chorionic Gonadotropin (β-HCG)	Myeloperoxidase (MPO)
B-Type Natriuretic Peptide (BNP)	Myoglobin
Creatine Kinase-MB (CK-MB)	N-Terminal Prohormone of Brain
D-Dimer	Natriuretic Peptide (NT-proBNP)
Digitoxin	Procalcitonin (PCT) ^(a)
High Sensitivity C-Reactive Protein (hsCRP)	Total Creatine Kinase (CK)
Homocysteine*	Troponin I
Human Chorionic Gonadotropin (HCG)	Troponin T

*Refer to Limitations of Procedure

STORAGE AND STABILITY

Unopened and opened vials of Omni•CARDIO stored at 2-8°C are stable as follows:

β-HCG:	15 days for Level UL, L, 1, 2, and 3
BNP:	15 days for Level UL, L, 1, 2, and 3
CK:	15 days for Level UL, L, 1, 2, and 3
CK-MB:	15 days for Level UL, L, 1, 2, and 3
D-Dimer:	15 days for Level UL, L, 1, 2, and 3
Digitoxin:	15 days for Level UL, L, 1, 2, and 3
HCG:	15 days for Level UL, L, 1, 2, and 3
hsCRP:	15 days for Level UL, L, 1, 2, and 3
MPO:	15 days for Level UL, L, 1, 2, and 3
Myoglobin:	10 days for Level UL, L, 1, 2, and 3
NT-proBNP:	10 days for Level UL, L, 1, 2, and 3
PCT:	15 days for Level UL, L, 1, 2, and 3
Troponin I*:	5-15 days for Level UL, L, 1, 2, and 3
Troponin T:	15 days for Level UL, L, 1, 2, and 3

Omni•CARDIO maintained at -25 to -15°C is stable until the expiration date as stated on the box.

* Refer to page 32 for lot-specific values of Troponin I's open-vial stability claim

SELF-DEFROSTING FREEZERS ARE NOT SUITABLE.

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

CONTROL RANGES

The published control ranges are based upon a combination of replicate assays of representative samples by participating laboratories, instrument/reagent manufacturers and direct correlation with other analytical systems in accordance with established protocol. Instrument values provided are specific to this lot of control only and are intended to assist the laboratory in establishing its own means and ranges. All values have been assigned with instruments and reagents available at the time of assay and expected values may vary with different reagents and/or methodologies. Laboratory established means should fall within the assigned ranges although subsequent instrument, reagent or calibration modifications may invalidate assigned values.

Peer comparison data and latest QC lot specific updates are available online through LabLink® xL Quality Assurance Program at www.maslablink.com. Refer to the Technical Assistance section for contact information.

INSTRUCTIONS FOR USE

Thaw control at room temperature (18-25°C) on a rocker or with periodic gentle inversion until liquid and then immediately store at 2-8°C. Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Open the vial and transfer the required quantity of control into a clean sample cup. Replace cap immediately and store the opened vial at 2-8°C. Assay controls in accordance with the reagent manufacturer's instructions for unknown specimens.

ONCE THAWED, DO NOT REFREEZE THE CONTROL.

Dropper tips are available to help minimize exposure of control samples to air. DO NOT use a syringe needle to withdraw sample through the control bottle cap.

QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURE

Compatibility of Omni•CARDIO has been demonstrated only with methods shown in this insert. Caution should be employed when using these controls with methods for which values have not been printed.

If low values are observed for Procalcitonin assay results upon initial thaw, an additional equilibration at 18-25°C for 4-6 hours may be necessary for full recovery.

Homocysteine is weighed into Omni•CARDIO; however, no claim is made for expected values nor the stability of this constituent. Testing of and assigning values to the Homocysteine constituent of this control is the responsibility of the individual laboratory.

If any constituent recovery falls below published ranges, contact MAS Technical Assistance for current peer group data. Peer comparison data and latest QC lot specific updates are available online through LabLink® xL Quality Assurance Program at www.maslablink.com.

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. This product is intended for use as an assayed control for quantitative assays of listed constituents in human serum. This product is not intended for use as a calibrator. For professional use only.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit www.maslablink.com and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit www.maslablink.com select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
OCRD-UL	MAS® Omni•CARDIO™, Level Ultra Low	6 x 3 mL
OCRD-L	MAS® Omni•CARDIO™, Level Low	6 x 3 mL
OCRD-101	MAS® Omni•CARDIO™, Level 1	6 x 3 mL
OCRD-202	MAS® Omni•CARDIO™, Level 2	6 x 3 mL
OCRD-303	MAS® Omni•CARDIO™, Level 3	6 x 3 mL
OCRD-MP	MAS® Omni•CARDIO™, Multi-Pack	6 x 3 mL
286-606	Dropper Tips	Pkg. 100

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2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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IVD

		CON									
		OCRD2605UA		OCRD2605LA		OCRD26051A		OCRD26052A		OCRD26053A	
			<R>		<R>		<R>		<R>		<R>
BECKMAN COULTER DXC / SYNCHRON SYSTEMS											
C-Reactive Protein, High Sensitivity (hsCRP)											
	Immunoturbidimetric (CRPH)	**		**		**		**		**	mg/dL
Creatine Kinase (CK)	Enzymatic	**		**		**		**		**	U/L
Creatine Kinase (CK)	Enzymatic (CK-)	**		**		**		**		**	U/L
Digitoxin	Latex Inhibition	**		**		**		**		**	ng/mL

		SI									
		OCRD2605UA		OCRD2605LA		OCRD26051A		OCRD26052A		OCRD26053A	
			<R>		<R>		<R>		<R>		<R>
BECKMAN COULTER DXC / SYNCHRON SYSTEMS											
C-Reactive Protein, High Sensitivity (hsCRP)											
	Immunoturbidimetric (CRPH)	**		**		**		**		**	mg/L
Creatine Kinase (CK)	Enzymatic	**		**		**		**		**	μkat/L
Creatine Kinase (CK)	Enzymatic (CK-)	**		**		**		**		**	μkat/L
Digitoxin	Latex Inhibition	**		**		**		**		**	nmo/L

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


LOT OCRD2605U
OCRD2605L
OCRD26051
OCRD26052
OCRD26053

2026-05-31

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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For insert updates go to:
www.thermoFisher.com/diagnostics

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
LOT OCRD2605U
OCRD2605L
OCRD26051
OCRD26052
OCRD26053

 2026-05-31

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
 scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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**** Data not available. If interested in participating in our value assignment process, please fax or email your contact information to our Value Assignment group at 510-771-1539, or mgc-va@thermofisher.com.**

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