# MAS® Omni•IMMUNE™

## LIQUID ASSAYED INTEGRATED IMMUNOASSAY CONTROL

**IVD**

### INTENDED USE

**Thermo Scientific MAS® Omni•IMMUNE™** is intended for use as an assayed control for monitoring assay conditions in many clinical laboratory determinations. Include **Omni•IMMUNE** with patient serum specimens when assaying for any of the listed constituents. Assay values are provided for the specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and instrument.

### PRODUCT DESCRIPTION

**Omni•IMMUNE** is a liquid stable control material prepared from human serum. Analyte levels are adjusted with various pure chemicals and preparations from human tissue or body fluids. Prostate specific antigen is from human seminal fluid, alpha-fetoprotein is from human cord serum and carcinoembryonic antigen is from hepatic metastases of human colon adenocarcinoma. Preservatives and stabilizers are added to maintain product integrity.

**CAUTION: Omni•IMMUNE** is prepared from human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found non-reactive for Hepatitis B Surface antigen (HBsAg), Hepatitis C (HCV), HIV-1 and HIV-2. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories, 2009." The packaging of this product contains dry natural rubber.

**DANGER: Omni•IMMUNE** control contains ≤0.1% bovine serum albumin (BSA) and ≤0.01% Proclin 300.

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

### CONSTITUENT LIST

17-alpha-OH-progesterone	Leuteinizing Hormone (LH)
25-Hydroxy Vitamin D	Lidocaine
Acetaminophen	Lithium
Adrenocorticotrophic Hormone (ACTH)	n-Acetylprocainamide hydrochloride (NAPA)
Alpha-Fetoprotein (AFP)	Prostatic acid phosphatase (PAP)
Aldosterone	PAPP-A*
Amikacin	Phenobarbital
Benzodiazepine (Nordiazepam)	Phenytoin
Beta - Human Chorionic Gonadotropin (β-hCG)	Phenytol, Free
Beta-2-Microglobulin	Primidone
Calcitonin*	Procainamide
Cancer Antigen 125 (CA 125)	Procalcitonin
Cancer Antigen 15-3 (CA 15-3)	Progesterone
Cancer Antigen 19-9 (CA 19-9)	Prolactin
Carbamazepine	Prostate Specific Antigen (PSA)
Carcinoembryonic Antigen (CEA)	Prostate Specific Antigen, free (Free PSA)
Cortisol	Parathyroid Hormone, Intact (PTH, Intact)
C-Peptide	Quinidine
Dehydroepiandrosterone- Sulfate (DHEA-Sulfate)	Salicylate
Digoxin	Triiodothyronine (T3)
Disopyramide	Triiodothyronine, free (Free T3)
Estradiol	Thyroxine (T4)
Estril, free	Thyroxine, free (Free T4)
Ethosuximide	Thyroxine binding globulin (TBG)
Ferritin	Testosterone
Folate	Theophylline
Fructosamine	Thyroglobulin
Follicle Stimulating Hormone (FSH)	Tobramycin
Gastrin	Tricyclic Antidepressants (TCA) (Nortryptiline)
Gentamicin	Thyroid Stimulating Hormone (TSH)
Growth Hormone	Thyroid Hormone Uptake (T-uptake)
Human Chorionic Gonadotropin (hCG)	Valproic Acid
Homocysteine	Valproic Acid, Free
Immunoglobulin E (IgE)	Vancomycin
Insulin*	Vitamin B12
Insulin-like Growth Factor - 1 (IGF-1)	

\*Refer to Limitations of Procedures

### STORAGE AND STABILITY

Unopened vials of **Omni•IMMUNE** are stable for 30 days from receipt when stored at 2-8°C. Once opened, vials of **Omni•IMMUNE** are stable for 30 days when stored tightly capped at 2-8°C. This product is stable until the expiration date on the box when stored at -25 to -15°C. **Self-defrosting freezers are not suitable.** If any constituent recovery falls outside the published ranges, contact MAS Technical Assistance for current peer group data. Peer comparison data and latest QC lot specific updates are available online through Lab Link xL Quality Assurance Program at [www.maslablink.com](http://www.maslablink.com).

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

### CONTROL RANGES

The published control ranges are based upon a combination of replicate assays of representative samples by participating laboratories, instrument/reagent manufacturers and direct correlation with other analytical systems in accordance with established protocol. Instrument values provided are specific to this lot of control only and are intended to assist the laboratory in establishing its own means and ranges. All values have been assigned with instruments and reagents available at the time of assay and expected values may vary with different reagents and/or methodologies. Laboratory established means should fall within the assigned ranges although subsequent instrument, reagent or calibration modifications may invalidate assigned values.

Peer comparison data and latest QC lot specific updates are available online through LabLink® xL Quality Assurance Program at [www.maslablink.com](http://www.maslablink.com). Refer to the Technical Assistance section for contact information.

### INSTRUCTIONS FOR USE

Thaw control at room temperature (18-25°C) on a rocker or with periodic gentle inversion until liquid and then immediately store at 2-8°C. Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Once control is removed from 2-8°C use immediately. Open the vial and transfer the required quantity of control into a clean sample cup. Replace cap immediately and store the opened vial at 2-8°C.

### ONCE THAWED, DO NOT REFREEZE THE CONTROL

Dropper tips are available to help minimize exposure of control samples to air. DO NOT use a syringe needle to withdraw sample through the control bottle cap.

### QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

### LIMITATIONS OF PROCEDURES

Compatibility of **Omni•IMMUNE** has been demonstrated only with methods shown in this insert. Caution should be employed when using these controls with methods for which values have not been printed.

Insulin, calcitonin and PAPP-A are weighed into **Omni•IMMUNE**; however, no claim is made for expected values nor the stability of these constituents. Testing of and assigning values to the Insulin and PAPP-A constituents of this control is the responsibility of the individual laboratory.

The specificity of antibodies used in immunoassay procedures may vary between lots. As with patient specimens, cross-reactivities between certain constituents in this control serum may occur with some reagents. Consult reagent manufacturer's package insert regarding possible interferences.

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. This product is intended for use as an assayed control for quantitative assays of listed constituents in human serum. This product is not intended for use as a calibrator. For professional use only.

### TECHNICAL ASSISTANCE

**In the USA**, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit [www.maslablink.com](http://www.maslablink.com) and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

**Outside of the USA**, if your laboratory subscribes to LabLink xL, visit [www.maslablink.com](http://www.maslablink.com) select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
OIM-101	MAS® Omni•IMMUNE™, Level 1	6 x 5 mL
OIM-202	MAS® Omni•IMMUNE™, Level 2	6 x 5 mL
OIM-303	MAS® Omni•IMMUNE™, Level 3	6 x 5 mL
286-606	Dropper Tips	Pkg. 100

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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
**LOT** OIM27031A  
OIM27032A  
OIM27033A

 2027-03-31

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

**MAS<sup>®</sup> Omni•IMMUNE<sup>™</sup>**


**LIQUID ASSAYED INTEGRATED IMMUNOASSAY CONTROL**

**IVD**

<div>U</div>	CON						SI							
	OIM27031A		OIM27032A		OIM27033A		OIM27031A		OIM27032A		OIM27033A			
	<div><math>\bar{x}</math></div>	<div>&lt;R&gt;</div>	<div><math>\bar{x}</math></div>	<div>&lt;R&gt;</div>	<div><math>\bar{x}</math></div>	<div>&lt;R&gt;</div>	<div>U</div>	<div><math>\bar{x}</math></div>	<div>&lt;R&gt;</div>	<div><math>\bar{x}</math></div>	<div>&lt;R&gt;</div>	<div><math>\bar{x}</math></div>	<div>&lt;R&gt;</div>	<div>U</div>
BIOMERIEUX VIDAS SYSTEMS														
Human Chorionic Ganadotropin (hCG) ELFA	**		**		**		mIU/mL	**		**		**		IU/L
Procalcitonin ELFA	**		**		**		ng/mL	**		**		**		µg/L
Prostate Specific Antigen (PSA) ELFA	**		**		**		ng/mL	**		**		**		µg/L
Prostate Specific Antigen, Free (fPSA) ELFA	**		**		**		ng/mL	**		**		**		µg/L
Thyroid Stimulating Hormone (TSH) ELFA	**		**		**		µIU/mL	**		**		**		mIU/L
Thyroxine, Free (FT4) ELFA	**		**		**		ng/dL	**		**		**		pmol/L

**LOT** OIM27031A  
OIM27032A  
OIM27033A

 2027-03-31

For Insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

**MAS<sup>®</sup> Omni•IMMUNE<sup>™</sup>**


**LIQUID ASSAYED INTEGRATED IMMUNOASSAY CONTROL**

**IVD**

CON						SI					
<div><div>U</div></div>											
OIM27031A		OIM27032A		OIM27033A		OIM27031A		OIM27032A		OIM27033A	
<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>	<div><div><div><div></div></div></div></div>
<div>U</div>						<div>U</div>					
BRAHMS KRYPTOR SYSTEMS											
Procalcitonin (PCT)											
TRACE-Technology											
**		**		**		**		**		**	
ng/mL						µg/L					

**LOT** OIM27031A  
OIM27032A  
OIM27033A

 2027-03-31

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

MAS® Omni•IMMUNE™

LIQUID ASSAYED INTEGRATED IMMUNOASSAY CONTROL

IVD

<div><div>U</div><div>CON</div></div>										<div><div>SI</div></div>									
		OIM27031A		OIM27032A		OIM27033A				OIM27031A		OIM27032A		OIM27033A					
		<div>X̄</div>	<div>&lt;R&gt;</div>	<div>X̄</div>	<div>&lt;R&gt;</div>	<div>X̄</div>	<div>&lt;R&gt;</div>	<div>U</div>		<div>X̄</div>	<div>&lt;R&gt;</div>	<div>X̄</div>	<div>&lt;R&gt;</div>	<div>X̄</div>	<div>&lt;R&gt;</div>	<div>U</div>			
DIASORIN LIAISON																			
C-Peptide	CLIA	**		**		**		ng/mL		**		**		**			nmol/L		
25-Hydroxy Vitamin D	Chemiluminescence	**		**		**		ng/mL		**		**		**			nmol/L		

**LOT** OIM27031A  
OIM27032A  
OIM27033A

 2027-03-31

**thermo**  
scientific







For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermoFisher.com/diagnostics](http://www.thermoFisher.com/diagnostics)

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 2027-03-31

**thermo**  
scientific

For insert updates go to:  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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

**If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.**

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**Microgenics Corporation**
**B-R-A-H-M-S GmbH**