MAS® T-Marker

LIQUID ASSAYED TUMOR MARKER CONTROL

IVD

INTENDED USE

MAS® T-Marker is intended for use in the clinical laboratory as a quantitative control for monitoring test procedures used to assay human serum for tumor markers. Include **T-Marker** with patient serum specimens when assaying for any of the listed constituents. Assay values are provided for the specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and instrument.

PRODUCT DESCRIPTION

T-Marker is a liquid stable control material prepared from human serum. Analyte levels are adjusted with various pure chemicals and preparations from human tissue or body fluids. Preservatives and stabilizers are added to maintain product integrity.

CAUTION: T-Marker is prepared from human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C (HCV), HIV-1 and HIV-2. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 2009. The packaging of this product contains dry natural rubber.

DANGER: T-Marker contains ≤0.2% bovine serum albumin (BSA).

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/ face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

STORAGE AND STABILITY

Unopened vials of **T-Marker** are stable for 90 days from receipt when stored at 2-8°C. Once opened, vials of **T-Marker** are stable for 30 days when stored tightly capped at 2-8°C. This product is stable until the expiration date on the box when stored at -25 to -15°C. **Self-defrosting freezers are not suitable.**

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

CONTROL RANGES

The published control ranges are based upon a combination of replicate assays of representative samples by participating laboratories, instrument/reagent manufacturers and direct correlation with other analytical systems in accordance with established protocol. Instrument values provided are specific to this lot of control only and are intended to assist the laboratory in establishing its own means and ranges. All values have been assigned with instruments and reagents available at the time of assay and expected values may vary with different reagents and/or methodologies. Laboratory established means should fall within the assigned ranges although subsequent instrument, reagent or calibration modifications may invalidate assigned values.

Peer comparison data and latest QC lot specific updates are available online through LabLink® xL Quality Assurance Program at www.maslablink.com. Refer to the Technical Assistance section for contact information.

INSTRUCTIONS FOR USE

Thaw control at room temperature (18-25°C) on a rocker or with periodic gentle inversion until liquid and then immediately store at 2-8°C. Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Once control is removed from 2-8°C use immediately. Open the vial and transfer the required quantity of control into a clean sample cup. Replace cap immediately and store the opened vial at 2-8°C.

ONCE THAWED, DO NOT REFREEZE THE CONTROL.

Dropper tips are available to help minimize exposure of control samples to air. DO NOT use a syringe needle to withdraw sample through the control bottle cap.

QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURE

Tumor Markers included in the control are AFP, CA 125, CA 15-3, CA 19-9, CA 27-29, CA 72-4, Calcitonin, CEA, β-hCG, PAP, PSA and Thyroglobulin.

The following constituents are also weighed into **T-Marker**. However, no claim is made for expected values nor the stability of these constituents: NSE, ACTH and Cyfra 21-1.

Compatibility of **T-Marker** has been demonstrated only with methods shown in this insert. Caution should be employed when using these controls with methods for which values have not been printed.

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. This product is intended for use as an assayed control for quantitative assays of listed constituents in human serum. This product is not intended for use as a calibrator. For professional use only.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit www.maslablink.com and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit www.maslablink.com select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
TUM-101	T-Marker, Level 1	6 x 3 mL
TUM-202	T-Marker, Level 2	6 x 3 mL
TUM-303	T-Marker, Level 3	6 x 3 mL
286-606	Dropper Tips	Pkg. 100

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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
LOT TM26051
TM26052
TM26053

 2026-05-31

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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
LOT TM26051
TM26052
TM26053

 2026-05-31

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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
LOT TM26051
TM26052
TM26053

 2026-05-31

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific







For insert updates go to:
www.thermofisher.com/diagnostics

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 2026-05-31

thermo
scientific

For insert updates go to:
www.thermofisher.com/diagnostics

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

If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.

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USA
☎ 800-232-3342
sales.diagnostics.fmt@thermofisher.com

China
☎ +86 800-810-5118
cdx.cn.info@thermofisher.com

New Zealand
☎ 0800 933 966
auinfo@thermofisher.com

Switzerland & Austria
☎ +41 26 663 86 70
cdx.ch.info@thermofisher.com

Asia Pacific
☎ +61 1800 333 110
cdd.asia.info@thermofisher.com

France
☎ +33 1 40 86 65 20
cdx.fr.info@thermofisher.com

Nordic
☎ +47 2 325 0433
info.nordic.cdd@thermofisher.com

United Kingdom & Ireland
☎ +44 1442 868 940
cdx.uk.info@thermofisher.com

Australia
☎ +61 1800 333 110
auinfo@thermofisher.com

Germany
☎ +49 0800-40 40 771
cdx.de.info@thermofisher.com

South Africa
☎ +27117926790
support-za.idd@thermofisher.com

Canada
☎ 800-282-4075
info.cddcanada@thermofisher.com

Japan
☎ +81 (0)120-147-075
JPYOK-CDD.QC@thermofisher.com

Spain, Portugal & Italy
☎ +34 93589 8338
cdx.es.info@thermofisher.com

For countries not listed:
distributor.cdd@thermofisher.com

Microgenics Corporation
46500 Kato Road
Fremont, CA 94538 USA
ISO 13485 Certified Company
www.thermofisher.com/diagnostics

B.R.A.H.M.S GmbH
Neuendorfstrasse 25
16761 Hennigsdorf, Germany
Tel: +49 (0) 800 404 077 12
Fax: +49 (0) 800 404 077 13

TUM-INS-VA
Rev. 21 2020 11

thermo
scientific